*Tecnovita*  
by **BH**

Instrucciones de uso y garantía  
Operating instructions and guarantee  
Mode d'emploi et garantie  
Instruções de utilização e garantia  
Manuale d'uso e garanzia  
Bedienungsanleitung und Garantie

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**902 170 258**

Thank you for purchasing our YR31 e-belt.  
All of our products undergo rigorous testing to ensure maximum quality.  
To use the unit correctly, first read through the whole manual before using it.

## **INTRODUCTION**


EMS, electric muscular stimulation is a widely recognized method that has been used on the sports field for a long time. The e-belt is both easy and comfortable to use. This product is specially designed to help tone the abdominal muscles. It allows you to enjoy the exercise and increase muscle capacity by developing strength and improving tone effortlessly. Thanks to the design of the bt belt electrostimulator it can tone the abdominal muscles anytime, anywhere.

## **GENERAL SAFETY INSTRUCTIONS**

1. Before using the appliance make sure it is in good working order. If you detect any faults get in touch with your local dealer.
2. The unit must not be used by children or the disabled without suitable adult supervision.
3. Pregnant women must not use it.
4. Do not use during the menstrual period.
5. Do not use if you suffer from heart disease or other type of disease.
6. Do not use the appliance if you have cancer.
7. Do not use if you are prone to attacks, for example, epilepsy.
8. Do not use the appliance after surgical operations in which increased muscle tension may affect the recovery process
9. Do not use directly on open wounds or varicose veins.
10. Do not use if you have any implanted electrical device, such as a pacemaker.
11. Do not use if you have any metal implants.
12. Do not use near the genitals.
13. Do not use on broken or inflamed skin.
14. Do not use the appliance near your heart.
15. The device is only for abdominal exercise.
16. If you suffer from any serious illness, especially if you suspect that you may have high blood pressure or problems with blood clotting, consult your doctor before using the appliance.
17. Consult your doctor before using this product if you have any type of skin disease.
18. If you suffer from chronic pain of unknown origin, regardless of where on your body, consult your doctor before using the appliance.

19. People undergoing medical treatment should consult a doctor before using this appliance.
20. Diabetics should consult their doctor before using this appliance.
21. People with constant skin irritation resulting from prolonged stimulation by the same electrode, should consult a doctor.
22. Consult your doctor if you have any problems resulting from the stimulation treatment.
23. Stop using the appliance if you experience any discomfort during use.
24. Do not use after consuming alcohol or drugs.
25. If the appliance malfunctions or if you feel discomfort or pain, discontinue use immediately.
26. Do not use while sleeping, driving or while using some other device.
27. Keep the appliance away from heat sources and do not use close (<1 m) to medium or microwave wave devices (e.g. mobile phones) because they can cause power surges.
28. Do not use when there are high frequency surgical devices connected simultaneously.
29. Do not allow the electrodes to come into contact with metal objects, e.g. belts or collars, during operation.
30. Do not use near flammable, gaseous or explosive substances.
31. Do not use half an hour before and after meal.
32. Switch the appliance off before removing it to prevent unnecessary irritation.
33. Do not apply enriching creams beforehand, as these cause the electrodes to wear out faster or may cause unwanted current surges.
34. New users should use this product starting from low frequency to avoid the possible uncomfortable feeling caused by stimulating muscles.
35. The maximum use time is 30 minutes.
36. Do not use the unit in damp/humid surroundings and prevent it from getting wet. Do not touch the device with wet hands while it is operating. Under no circumstances must this product be used in the bath or shower.
37. Do not expose the appliance to direct sunlight or high temperatures.
38. Do not allow the appliance to come into contact with water or any other liquid. Do not immerse the appliance in liquid to clean it.
39. If you notice that the appliance is not operating correctly, do not attempt to repair, dismantle or modify it. Unauthorised repairs may cause injury, in addition to breaking the terms of the warranty. Contact your technical assistance service.

## DESCRIPTION


- 1 – BELT
- 2 – EXTENSION BAND
- 3 – CONTROLLER
- 4 – BATTERIES AAA LR03
- 5 – SQUARE SILICONE PATCH
- 6 – RECTANGLE SILICONE PATCH X 2


## OPERATING INSTRUCTIONS

### BEFORE USE

1. To insert batteries correctly, open the back battery cover and insert 2 alkaline batteries AAA in the battery holder according to the positive and negative terminals. Close the battery holder.

- Insert the small head of the connecting line into the socket on the top of the controller. Then insert the larger head of the connecting line into the connecting socket of waistband,


- Remove the protective film. Note that it is needed to remove the protective film of the side printed with yellow gridding.


- Paste the silicone patches at the right position and refer to the positions of the patches at the sides:


M-L	Waist 68-82 cm
X-L	Waist 84-102 cm
XX-L	Waist 104-130 cm

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
6. Fix the belt to your body, the silicone patches should be stuck on the skin. Use extension band if necessary.

## REMOTE CONTROL

- 1 – ON/OFF button – Press for turn ON/OFF
- 2 – Intensity +
- 3 – Intensity -
- 4 – Select mode (1-8)
- 5 – Intensity (1-18)
- 6 – Mode
- 7 – Time

NOTE: The recommended time of use is 30 minutes per day. If you wish to use it again, wait 60 minutes and then turn it on.

## **MODES**

Mode 1: Normal sit-up mode:

Exercise anterior rectus abdominis and transverses abdominis to strengthen the core of abdomen and share the work of muscle group of waist, so as to protect lumbar vertebrae, improve body balance, increase basal metabolic rate, promote heat consumption, reduce abdominal fat and reach the effect of body shaping.

Mode 2: Normal abdominal curl with trunk twisting mode:

Exercise intra-abdominal external oblique muscles at the both sides of abdomen, strengthen gravity center control of body rotation and lateral bending and improve the stability of body.

Mode 3: Abdomen extension mode:

Make centrifugal movement of muscles of upper and lower abdomen to stretch and relax abdominal muscles to reach the shaping effect of abdominal muscles.

Mode 4: Sit-up speed mode:

Exercise the whole abdominal muscle group (surface muscle and deep muscle) through fast-slow inputting of sports signal, so as to reach the effect of aerobic exercise of abdomen, meanwhile, promote core stability and strengthen body control and balance.

Mode 5: flank bending and extension mode:

Make abdomen have the effect of aerobic exercise and improve the control capacity of flank through orderly exchange movement of flank muscles

Mode 6: Abdominal core strengthening:

Strengthen to exercise abdominal deep muscle group (internal oblique muscle and transverse muscle of abdomen), enhance abdominal core control and improve the balance force of body.

Mode 7: Abdominal extension and relaxation:

Relaxation massage of the whole abdominal muscle group can relieve the fatigue of muscles.

Mode 8: Synthetic abdominal sports:

This mode will exercise the whole abdominal surface and deep muscles, strengthen core control and protect back muscle group and vertebra. Besides, it enables body bending forward, twisting and lateral bending much more flexible, to improve the abdominal muscle and to achieve abdomen tightening and WC reduction.

## USE REFERENCE

Week	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
1	Intensity 1	Rest	Intensity 2	Rest	Intensity 4	Intensity 4	Intensity 4
2	Rest	Intensity 4	Intensity 4	Rest	Intensity 8	Intensity 8	Intensity 8
3	Rest	Intensity 9	Intensity 9	Rest	Intensity 9	Intensity 10	Intensity 10
4	Rest	Intensity 9	Intensity 9	Rest	Intensity 9	Intensity 9	Intensity 10

This is a reference plan of use. For make your personal plan you should take into account your health and personal conditions.

## CLEANING AND MAINTENANCE

1. After using wipe gently the stains such as perspiration on belt with a slightly damp cloth.
2. The waistband could only be washed by hands. Remove the controller and silicone patch before washing the product, besides, the silicone patch shall be pasted with protective film.
3. Please plug the interface of connecting line on the waistband with stopper when cleaning.
4. Clean the waistband by hands with warm water, do not use washing machine.

5. Do not use whitening or cleaning products.

6. Leave the item to get dry on air. Do not use dryers. Do not iron the item.
7. Store the waistband at a dry place.

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1. Silicone patch shall be used by one person to keep healthy.
2. Please clean the skin that will contact this product before using.
3. Avoid the contact of the fingers with the patches.
4. The surface of silicone patch shall be clean to avoid dust, oily articles, sticky articles and other dirt, which may make the viscosity decrease; please use a little flowing water to clean and fingers to wipe the dirt on the surface of silicone patch, after that, hang to dry the silicone patch to recover its viscosity. However if much more water is used when washing, the viscosity will also decrease.
5. Please do not wash the silicone patch frequently, or its service life will be shortened.
6. Please paste the transparent plastic protective film again after use and put the silicone patch at the shady and cool place for storage.
7. Since the silicone patch belongs to consumables, it is advised to consider purchasing new one when its viscosity could not be recovered or the stimulus decreases after washing many times.

## TECHNICAL SPECIFICATIONS

MODEL	E BELT YR31
FREQUENCY	1Hz – 1000Hz
IMPULSE WIDTH	20µs - 400µs
BATTERIES	2X LR03 AAA ALKALINE BATTERIES
DIMENSIONS BELT	1040X172MM
EXTENDABLE BELT	300X70MM
WEIGHT	54g

### Warranty

This product is guaranteed for 2 years in accordance with that stated in Directive 1999/44/EC of the European Parliament and of the Council 25<sup>th</sup> May, 1999, regarding certain aspects of the sale of consumer goods and associated guarantees, incorporated into Spanish Law by virtue of Law 23/2003, 10<sup>th</sup> July; Guarantees on the sale of consumer goods.

To this effect, EXERCYCLE, S.L. VAT n° B-01029032, with business address c/ Zurrupitieta 22, Vitoria, Spain, guarantees this product.

<b>Date of purchase.-</b>	<b>Purchased in.-</b>
<b>Purchaser.-</b>	<b>Dealer stamp.-</b>

Please keep the receipt as proof of purchase.

Should you have any queries please do not hesitate to phone our customer service line:

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BH China Co., Ltd.